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The Kerbl quick test for the qualitative detection of rotaviruses, bovine corona viruses (BCV), E.coli-K99 and cryptosporidium parvum in the faeces of calves.

## Application

A) Pick up the stool with the pipette supplied.

### Explanatory Note:

- For solid stools use the applicator in the cap of the sample container to pick up samples. In doing so, ensure that sufficient stool is picked so that the test fluid becomes clearly discoloured.
- If possible, avoid picking up mucus.
- Stool samples from the stall floor can falsify the results.

B) Using the pipette, put 3-5 drops of the stool into the sample container. The test fluid must become clearly discoloured.


C) Close the sample container and shake it lightly.

D) Take the test cassette out of the sealed aluminium bag and place the cassette on a level surface.

E) Unscrew the cap from the tip of the sample container.

F) Put 3 drops of the test mixture in each of the test cassette sample well.


### Note:

- If the fluid does not raise the test strip, additional drops may be put in the sample well.
- No fluid must drip onto the display panels.


G) After 10 minutes read the result in the display panels.

A **positive test result** is shown when the T and C markings are visible.

A **negative test result** is shown when only the C marking is visible.

### Storage

The test must be stored and carried out in a dry location and at a temperature between 2 °C and 30 °C.

### Liability

The complete liability risks in relation to the use of this product are borne by the purchaser. The manufacturer assumes no liability for indirect, special or consequential damages of any kind whatsoever which result from the use of, or the carrying out and evaluation of tests with this product.

### Disposal

No special means of disposal are necessary.

### Package Contents

- 5 x Test cassettes in shrink-wrapped aluminium bags
- 5 x Sample containers with test fluid
- 5 x Pipettes
- 1 x Instructions

### Please Note:

Test for single use.  
After opening the test cassette bag carry out the test as soon as possible. The test cassette packing must not be damaged prior to use. Observe the evaluation time.

### Sensitivity and Specificity of the 4D-Test:

	Sensitivity	Specificity
Rotaviruses	99,0%	98,0%
Bovine coronaviruses	98,4%	98,0%
E.coli-K99	97,8%	99,0%
Cryptosporidium parvum	98,2%	99,0%
Study vs. ELISA		

Albert Kerbl GmbH  
Felzenzell 9  
84428 Buchbach  
Germany  
Tel. +49 8086 933-100  
Fax +49 8086 933-500  
info@kerbl.com  
www.kerbl.com

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